Sattelprothese Endo-Modell®
Endo-Model® Saddle Prosthesis

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CE 0123

Waldemar Link GmbH & Co. KG

Barkhausenweg 10 · 22339 Hamburg, Germany

P.O. Box 63 05 52 · 22315 Hamburg, Germany

Tel.: +49 40 53995-0 · Fax: +49 40 5386929

E-mail: info@linkhh.de · Internet: www.linkhh.de

Sattelprothese Endo-Modell®

Endo-Model® Saddle Prosthesis

Systembeschreibung	System Description
02 Problemstellung	The Problem
03 Systemübersicht	System Overview
 05 Indikationen/Kontraindikationen	 Indications/Contraindications
 Implantate	 Implants
06 Sattelprothese	Saddle Prosthesis
08 Prothesenschaft SPII® Modell Lubinus®	Lubinus SPII® Prosthesis Stem
 10 Instrumente	 Instruments
 OP-Technik	 Surgical Technique
12 Allgemeines	General
14 Spezielles: Dysplasietyp	Specifics: Dysplasia Type
15 Spezielles: Protrusionstyp	Specifics: Protrusion Type
16 Spezielles: Resektionstyp	Specifics: Resection Type
16 Postoperative Behandlung	Postoperative Management
17 Fallbeispiele	Case Histories
 Zubehör	 Accessories
20 Sonderanfertigungen	Custom-made Implants
20 Röntgenschablonen	X-ray Templates
 21 Literatur	 Literature
 24 Index	 Index
 Wichtige Hinweise	 Important Information

■ Problemstellung

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■ The Problem

In patients with extensively damaged or resected pelvic structures, often only the superior portions of the ilium are preserved, with remnants of the anterior or posterior cortical bone. Occasionally, there may even be loss of continuity with the lower parts of the pelvis. Since the superior portions of the ilium extend as far as the sacroiliac joint and the cortical bone thickens at the apex of the flared pyramidal remnant above the acetabulum, the surface merges in the superomedial region. More laterally, the cortical bone gets progressively thinner. Consequently, only the thick pyramid-shaped bone above the acetabulum and the thicker medial area of the iliac wing can be used for implant bearing.

The situation described above generally arises in cases of revision arthroplasty or tumour surgery where the acetabulum can not be reconstructed, either because a graft can not be placed or obtained or because the patient's condition rules out such a time-consuming operation.

Such cases nearly always require comprehensive surgical planning, for which computed tomography, magnetic resonance imaging, intravenous pyelography and 3D modelling (LINK® ComForm) all offer insights.

The LINK® Endo-Model® Saddle Prosthesis is an effective alternative in such cases.

The Saddle Prosthesis was developed in 1979 by Nieder et al. known as the Endo-Model® Mark I, it was a one-piece femoral stem with a crutch-shaped saddle. The model underwent further development in 1987, resulting in a saddle prosthesis with axial rotation known as the Endo-Model® Mark II.

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Abb./Fig. A

■ System Overview

The LINK® Endo-Model® Mark II Saddle Prosthesis consists of a symmetrical, saddle-shaped endpiece which articulates directly with a notch in the ilium. This saddle straddles the notch in the bone, the medial horn lying inside the pelvis and the lateral horn outside. The rotatable saddle component is seated on the post of the respective base component. Located between the base component and the saddle component is a hat-shaped UHMWPE sleeve which allows low-friction movement between these two chrome-cobalt components. A screw, which passes through the butt of the saddle to engage a circular channel running around the post of the base component, with-out touching the post metal, helps to prevent joint dislocation (**Fig. A**).

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The base component has a female taper distally, enabling it and the saddle (instead of the customary ball head) to be mounted on the chosen femoral stem component. The chosen stem components must be designed to fit the base components. Such stem components have a 14/16 mm taper, the interface of which has an eccentric hole that engages the pin in the base of the female taper of the base component, thus preventing rotation. The base component is additionally secured to the femoral component taper by means of at least two fixation screws which are screwed into two of the four distal threaded holes in the base component and into the predrilled recesses on the taper of the femoral component.

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The articulation between saddle and base component allows all the requisite rotational movements not permitted by the saddle joint itself. The pivoting movements that occur in the joint between bone and saddle as well as the rotation between the saddle and the base component replace the movements permitted by the “lost” ball-and-socket joint (**Fig. B**).

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The base components are angled in the shape of a bayonet and displace the femur laterad. This arrangement also permits standing on one leg, providing sufficient abductors are still present (**Fig. C**).

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
■ Indications/ Contraindications

Note: For specific indications/
contraindications, see page 23.

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
The Saddle Prosthesis may only be combined with LINK® femoral stems featuring a 14/16 mm taper and a hole at the interface to engage the anti-rotation pin inside the base component.

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The base component must be additionally secured to the taper with the supplied set screws. The appropriate tools are provided (**F** to **H** listed below). To achieve a secure screw fixation, the twist drill and the drill sleeve must be used to prepare blind holes (min. 2, max. 4) into the prosthesis taper.

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
For custom-made saddle component, see page 20.

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Saddle Prosthesis, complete, incl. 4 set screws

Art.-Nr. Item No.	Größe Size	Länge Length (mm)
15-8007/01	I	60
15-8007/02	II	70
15-8007/03	III	80
15-8007/04	IV	90
15-8007/05	V	100
15-8007/06	VI	130
15-8007/07	VII	160

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* Base components without anti-rotation pin for femoral stems with 12/14 mm taper available on request.

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

The Endo-Model® Saddle Prosthesis may be combined with 14/16 mm taper „Lubinus SPII® Prosthesis System“ femoral stems with holes for antirotation pins. Information on possible combinations with other LINK® femoral stems available on request.

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Prosthesis Stems with Bore Hole 14/16 mm taper for Saddle Protheses, Lubinus SPII®

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The size of the femoral stem, composed of a letter and a number, is indicated on the top of the collar:

R = rechts/right	3 = stark/large
L = links/left	4 = extrastark/x-large
01 = extraschlank/x-narrow	45 = extrastark plus/ x-large plus
1 = schlank/narrow	5 = superstark/super-large
2 = normal/medium	

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
To displace the intrapelvic soft tissues craniad, ventrad and dorsad.

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Ruskin rongeur, modified for gradual opening of the acetabular floor to take the medial horn of the saddle, 285 mm

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The lateral acetabular rim of the ilium is notched for the saddle using the bone rongeur and/or a long osteotome.

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The Lambotte osteotomes **15-2574/04** and **15-1574/06** and the gouge **15-1074** are used to open the floor of the acetabulum step by step. The inner and outer cortical bone at the centre of the original acetabular floor is carefully perforated to accommodate the medial saddle horn.

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Impactor for saddle base component, 220 mm

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Screwdriver for fixation screw (saddle) and set screw (base component), 210 mm

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
The holding forceps for holding the set screws have a female thread at the tip. This enables the surgeon to position and screw in the screw while holding it with the holding forceps.

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■ General

The different types of bone loss in the pelvic region can be broadly classified as follows:

1. **Dysplasia**
2. **Protrusion**
3. **Resection**

In dysplasia, bone loss is generally due to cranio-lateral migration or loosening of an artificial acetabular component. In the protrusion type, bone loss may also be caused by cranio-medial migration. It should be noted that protrusion may be either separated from the inner pelvis or open. Bone defects characterised by an open protrusion, with an interruption of the anterior and posterior connection, are usually the result of tumour resection (**Fig. D**).

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
Generally speaking, the chosen base component length should facilitate reduction by working with the remaining soft-tissue structures, given sufficient tension. This type of soft-tissue balancing helps reduce the chance of dislocation.

In revision arthroplasty, short base components are used in most cases and soft-tissue scarring often poses a problem for reduction. In contrast, in bone-tumour surgery, the longer base components tend to be used, since the articulation level is generally higher and there is usually little tension in the normal muscles that remain after resection; leg lengthening may be unavoidable. If the patient has a high risk of dislocation – this evaluation may be done intraoperatively – it may be necessary to limit critical movements for a few weeks. If necessary, the patient is immobilised by applying either light traction or a Whitman's plaster after the wound has healed.

The risk of dislocation becomes minimal once a fibrous pseudocapsule has formed around the saddle joint after a few weeks, which then continues to ossify.

If, however, dislocation does occur, an open reduction is nearly always required, since the saddle twists when a closed reduction is attempted and the medial horn makes contact with the bone.

When attaching the saddle component to the base component, it is important to ensure that the fixation screw, which fixes the saddle to the base component, is inserted from lateral so that it is accessible in the event of revision.

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
■ Specifics: Dysplasia Type

In dysplasia, the cranio lateral acetabular rim is no longer present. The floor of the acetabulum is usually closed and, therefore, needs to be opened in order to introduce the saddle. The triangular retrograde angled raspatory is then used to displace the intrapelvic soft tissues cranial, ventrad and dorsad. The acetabulum must then be resected so that any dorsal bone structures which could impede saddle pivoting are removed. In the dysplasia type, the area of saddle articulation tends to be medial.

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
■ Specifics: Protrusion Type

In the case of protrusion, the lateral acetabular rim, which is generally preserved, is notched to form a bearing surface for saddle articulation. If the intra-pelvic area is separated by a (generally thin) lamella of bone, the latter must be perforated to accommodate the medial horn of the saddle within the pelvis. It is important to ensure that the saddle horn is not located between the inner and outer cortical bone. In the case of protrusions without such separation, only the acetabular rim is notched, if necessary. Any caudal bone portions that could impede saddle mobility are removed in this case too. In the protrusion type, the site of saddle articulation is generally more lateral.

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■ Specifics: Resection Type

In the resection type, the inferior margin of the ilium is notched to the approximate depth of the saddle centre. The bone nibbling rongeur **15-1075** and the extra-long Dahmen gouge **75-3610/20** are ideal for notching the ilium.

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■ Postoperative management

After-treatment is similar to that for total hip prostheses. The patient does isometric muscle contraction, foot rotation and breathing exercises on the first postoperative day. Once the wound drainage devices have been removed and the dressings changed, the patient may stand up with the assistance of a physical therapist and take his first steps with the aid of a walker, at which point he may put weight on the operated leg. Mobilisation proceeds with physio-therapy exercises with the aid of a walker, elbow crutches or cane. Plenty of support should be provided – long-term if the patient shows unsteadiness. In tumour cases involving extensive resection in particular, it must be determined whether restriction of movement or even immobilisation is necessary for a time.

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■ Case Histories

Dysplasia type

A 66-year-old female patient who, due to dysplastic coxarthrosis, received a total hip replacement on the right in 1969 and on the left in 1972.

Revision surgery on the left was performed in 1988 due to loosening. In 1989, the left hip replacement showed major failure again.

The left hip replacement was revised in 1989; an exchange operation was performed. A saddle prosthesis was implanted due to severe destruction of the acetabular bone stock resulting from dysplasia. After the bone had been perforated, the head of the saddle was supported on the iliac wing. The saddle was connected to an Endo-Model® femoral stem size 5 via a no. 1 base component.

In 1990, there was evidence of good incipient sclerosis in the weight-bearing area of the iliac wing following slight cranial migration of the saddle head.

By 1995, the progression of sclerosis in the weight-bearing area of the saddle head was such that the medial horn was completely embedded in the bony bearing.

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Protrusion type

A 72-year-old male patient had received a right total hip replacement in 1975 due to coxarthrosis. A revision was performed in 1983 following periprosthetic infection. It was subsequently necessary to replace the acetabular component twice due to loosening.

The patient again showed loosening of the acetabular component, accompanied by destruction of acetabular bone stock and protrusion.

In 1988, revision surgery was performed on the right hip using a 200 mm-long revision stem, on which a medium-length base component with saddle was mounted.

In the protrusion type, the lateral acetabular rim, which is generally preserved, is notched for saddle articulation. If the inner pelvic area is separated by a thin lamella of bone, the latter must be perforated to accommodate the medial horn of the saddle within the pelvis. In the case of open protrusions, the acetabular rim is notched, if necessary.

Since 1992, the position of the saddle in relation to the pelvis has remained unchanged, with a good zone of sclerosis.

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
Resection type

A chondrosarcoma is tentatively diagnosed in a 46-year-old patient.

The preoperative X-ray suggests a tumour of the right acetabulum.
A chondrosarcoma is diagnosed on the basis of a biopsy.

Postoperative analysis of the specimen confirms the preoperative diagnosis.

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The postoperative X-ray shows the implantation of an acetabular component into a large homologous pelvic graft, which was fixed with screws and a plate.

The X-ray status three months after surgery shows that the reconstruction has collapsed.

Postoperative X-ray three months after implantation of a saddle prosthesis. Incipient reossification is discernible around the lateral horn of the saddle.

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
Saddle component with suture hole,
available on request

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
for LINK® Endo-Model® Saddle Prosthesis,
110% actual size, 2 sheets

E. Nieder, R. A. Elson, E. Engelbrecht, M. R. Kasselt, A. Keller, K. Steinbrink

The Saddle Prosthesis for salvage of the destroyed acetabulum
J Bone Joint Surg (Br) 1990, 72-B:1014-22 [H 67](#)

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Reconstruction using the saddle prosthesis following excision of primary and metastatic periacetabular tumors
Reprinted from Clinical Orthopaedics May 1995, Vol. 314 [H 77](#)

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The saddle prosthesis in pelvic primary and secondary musculoskeletal tumors:
functional results at several postoperative intervals
Springer Verlag, Berlin, 2000 [H 98](#)

Rubert, Malawer, Rosenthal

Reconstruction With the Saddle Prosthesis: Indications, Surgical Technique, and Results
Reprinted from Seminars In Arthroplasty, Vol. 10, Nr. 3 [H 100](#)

J. Benevenia, F. P. Cyran, J. S. Biermann, F. R. Patterson, M. C. Leeson

Treatment of Advanced Metastatic Lesions of the Acetabulum Using the Saddle Prosthesis
Clinical Orthopaedics & Related Research No. 426, pp. 23–31, 2004 [H 136](#)

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■ Indications/Contraindications

	Products		
	MP® Reconstruction System	MEGASYSTEM-C®	Endo-Model® Saddle Prosthesis
General Indications			
Mobility-limiting diseases, fractures or defects of the hip joint and the proximal femur which cannot be treated by conservative or osteosynthetic procedures	X		X
Mobility-limiting diseases, fractures or defects of the hip joint, the proximal and distal femur through the proximal tibia which cannot be treated by conservative or osteosynthetic procedures		X	
Indications			
Revision arthroplasty due to juxta-articular bone defects	X	X	X
Revision of loosened femoral prosthesis components involving extensive bone resorption of the proximal femur and widening of the medullary cavity or marked thinning of proximal femoral cortical bone	X	X	
Revision of loosened femoral prosthesis components by peri-/subprosthetic fracture	X	X	
Deformed proximal femur due to fractures or osteotomies	X	X	
Correction of bone deficiencies, e.g. due to tumours	X	X	X
Large post-revision and post-trauma segmental bone defects	X	X	X
Oncological and revision surgery from tibial to hip area (in conjunction with Endo-Model® SL® Rotational and Hinge Knee Prostheses)		X	
By severe or non-reconstructive acetabular defects			X
Contraindications			
Acute or chronic infections, local or systemic	X	X	X
Allergies due to (implant) materials	X	X	X
Revisions in septic environment	X	X	X
For preparation of the prosthesis bearing insufficient length of intact diaphysis (less than 80 mm)	X	X	
Distinctive muscular, nerve, vascular or other diseases which put the affected limb at risk	X	X	X
Insufficient bone integrity which prevents a stable anchorage of the prosthesis	X	X	X
Non preserved abductor- or iliopsoas musculature			X
Relative Contraindications			
Adiposity	X	X	X
Lacking of or foreseeably unassured compliance	X	X	X
Foreseeable overload/overstressing of joint prosthesis	X	X	X

Please note: These indications/contraindications refer to standard cases. The ultimate decision on whether or not an implant is suitable for a patient must be made by the surgeon based on his/her individual analysis and his/her experience.

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Please note the following regarding the use of our implants:

- 1. Choosing the right implant is very important.**
The size and shape of the human bone determine the size and shape of the implant and also limit the load capacity. Implants are not designed to withstand unlimited physical stress. Demands should not exceed normal functional loads.
- 2. Correct handling of the implant is very important.**
Under no circumstances should the shape of a finished implant be altered, as this shortens its life span. Our implants must not be combined with implants from other manufacturers.

The instruments indicated in the Surgical Technique must be used to ensure safe implantation of the components.
- 3. Implants must not be reused.**
Implants are supplied sterile and are intended for single use only. Used implants must not be reused.
- 4. After-treatment is also very important.**
The patient must be informed of the limitations of the implant. The load capacity of an implant cannot compare with that of healthy bone!
- 5. Unless otherwise indicated, implants are supplied in sterile packaging.**
Note the following conditions for storage of packaged implants:
 - Avoid extreme or sudden changes in temperature.
 - Sterile implants in their original, intact protective packaging may be stored in permanent buildings up until the "Use by" date indicated on the packaging.
 - They must not be exposed to frost, dampness or direct sunlight, or mechanical damage.
 - Implants may be stored in their original packaging for up to 5 years after the date of manufacture. The "Use by" date is indicated on the product label.
 - Do not use an implant if the packaging is damaged.
- 6. Traceability is important.**
Please use the documentation stickers provided to ensure traceability.
- 7. Further information** on the material composition is available on request from the manufacturer.

Follow the instructions for use!

Waldemar Link GmbH & Co. KG, Hamburg

All content in this catalog, including text, pictures and data, is protected by copyright. Every instance of use, whether in part or in whole and which is not permitted by the copyright act, is subject to our prior consent. In particular, this applies to the reproduction, editing, translation, publishing, saving, processing, or passing on of content stored in databases or other electronic media and systems, in any manner or form. The information in the catalogs is solely intended to describe the products and does not constitute a guarantee.

The Surgical Technique described has been written to the best of our knowledge and belief, but it does not relieve the surgeon of his/her responsibility to duly consider the particularities of each individual case.

Unless otherwise indicated, all instruments are made of surgical stainless steel.

Waldemar Link GmbH & Co. KG

Barkhausenweg 10 · 22339 Hamburg, Germany
P.O. Box 63 05 52 · 22315 Hamburg, Germany
Tel.: +49 40 53995-0 · Fax: +49 40 5386929
E-mail: info@linkhh.de · Internet: www.linkhh.de

LINK[®]